

At a glance

# Gas permeable and hybrid lenses

GAS PERMEABLE & HYBRID											
Minus BVP Range	PRODUCT NAME	BOZR mm	TD mm	H <sub>2</sub> O %	CT mm @ -3.00D	TINT	Dk	CYLINDERS	AXES	PRISM	Plus BVP Range

SPHERES & ASPHERES											
-30	POLYCON*	7.00 to 8.60	9.00 9.50		0.13' 0.12 <sup>2</sup>	Var.	12 to 151				+30
	ASPHERICON*	7.00 to 8.60	9.30 9.80		0.12	Var.	12 to 95				+30
-30	ASTROCON*	7.00 to 8.60	9.60 10.20		0.135	Var.	12 to 95				+30
-30	K DESIGN*	7.00 to 8.60	9.00 9.50		0.17 0.15	Var.	12 to 151				+30
-30	KE DESIGN*	7.00 to 8.60	9.20 9.80		0.15 0.16	Var.	12 to 151				+30
-30	PRACTITIONER SPECIFIED*	to Rx	to Rx		to Rx	Var.	12 to 151				+30
-25	AQUILA®	7.30 to 8.50 7.40 to 8.60	9.30 9.80 10.30		0.18	Light Blue	143				+25
	PERSECON® E	6.90 to 8.60 7.00 to 8.70 7.10 to 8.80 7.20 to 8.90	8.80 9.30 9.80 10.30		0.17	Clear or Grey	8				+25
-25	PERSECON® 92E	7.20 to 8.50 7.30 to 8.50 7.40 to 8.60	8.80 9.30 9.80 10.30		0.17	Light Blue	92				+25
	PERSECON® CE	7.20 to 8.40 7.30 to 8.50 7.40 to 8.60	8.80 9.30 9.80 10.30		0.17	Light Blue	54				+25
-30	CONTINUOUS CURVE (PMMA)	7.00 to 8.60	9.00 9.50		0.17 0.20	(912) Grey	0				+30

KERATOCONUS											
-30	SOPER KERATOCONUS*	to Rx	to Rx		to Rx	Var.	12 to 65				+25
-30	WOODWARD KERATOCONUS*	to Rx	to Rx		to Rx	Var.	12 to 65				+25
-25	PERSECON E KERATOCONUS	5.90 to 7.00 6.00 to 7.70	9.30 9.80		0.203	Clear or Grey	8				+25

Coding key for powers ■ 0.25 steps ■ 0.50 steps

GAS PERMEABLE & HYBRID											
Minus BVP Range	PRODUCT NAME	BOZR mm	TD mm	H <sub>2</sub> O %	CT mm @ -3.00D	TINT	Dk	CYLINDERS	AXES	PRISM	Plus BVP Range

TORIC											
-30	FRONT SURFACE TORIC	to Rx 7.00 to 8.60 (0.05)	to Rx 8.00 to 10.00		to Rx Various	Var.	12 to 65 Var.	up to 3.00D	1° ATC	up to 3.00Δ	+30
-30	BACK SURFACE TORIC	to Rx 7.00 to 8.60 (0.05)	to Rx 8.00 to 10.00		to Rx Various	Var.	12 to 65 Var.	up to 5.00D			+30
-30	BI-TORIC	to Rx 7.00 to 8.60 (0.05)	to Rx 8.00 to 10.00		to Rx Various	Var.	12 to 65 Var.	up to 5.00D			+30
-30	TORIC PERIPHERY	to Rx	to Rx 8.00 to 10.00		to Rx Various	Var.	12 to 65				+30

HYBRID												
-13	-8	SOFTPERM	7.10 to 8.10	14.30	25	0.11	Clear	14&550			+4	+6

MULTIFOCAL												
	-30	ASTROCON MULTIFOCAL	7.10 to 8.60	8.00 to 10.00		0.20	Blue or Grey	35 or 58	UP TO +3.50		+30	

Coding key for powers ■ 0.25 steps ■ 0.50 steps

For assistance with any contact lens fitting, please call Professional Services – 01489 775 533

**\*These designs are available in the following materials**

	Dk
Polycon II	12
Fluorocon 30	35
HDS	58
Fluorocon 60	65
Fluorocon 90	95
Fluorocon 151	151

1 Polycon lens design lenses made in Fluorocon 90 and Fluorocon 151 are 5/100th thicker (0.18 @ -3.00)

2 Polycon lens design lenses made in Fluorocon 90 and Fluorocon 151 are 5/100th thicker (0.17 @ -3.00)

See fitting guides on pages 37-38.

**Warranty conditions – made to order, traditional soft and RGP lenses**

All traditional soft and gas permeable lenses can be exchanged or returned for up to 90 days from the date of the **original order** for each patient, for the purpose of fitting, vision improvement or unsuccessful fits.

Only lenses returned within this time with the relevant order number will be credited to your account.